

A comprehensive method for the screening of drugs in severely decomposed human tissues by fast gas chromatography/tandem mass spectroscopy (Fast GC-MS/MS)

Brian Waters

Fukuoka University, Japan

Abstract

Statement of the Problem: Drug screening is an important reference in forensic autopsy investigations. In post-mortem toxicology, often the samples provided for analysis are in a severe state of putrefaction or decomposition. The presence of breakdown products such as lipids and amino acids make extraction of the compounds of interest difficult. Also, developing an analytical method capable of detecting trace levels of analytes from the interfering substances present in these complex matrices adds to the challenges.

Methodology & Theoretical Orientation: For this study, putrefied and decomposing tissue samples from actual cases autopsied at our department were analyzed. Human tissue specimens consisted of liver, kidney, spleen, lung, muscle and brain, if available. The drugs detected from these specimens included phenobarbital, chlorpromazine, promethazine, aripiprazole, amlodipine, telmisartan, rosvastatine, chlorpheniramine, etizolam and zolpidem. Specimens of 0.3 g were treated with urease, acidified or alkalized and extracted with acetonitrile. Lipid-removal and solid-phase extraction cartridges were employed while carefully monitoring the pH of samples to ensure the adequate removal of interfering substances. The extracts were evaporated and reconstituted in n-propyl acetate:methanol (1:1) for fast GC-MS/MS analysis.

Findings: The developed method was successful in clearly identifying drugs from putrefied specimens. The use of tandem mass spectrometry helped to reduce the influence of background noise and interfering substances.

Conclusion & Significance: Putrefied specimens are often the only remaining samples left from severely decomposed cadavers. The combination of a robust preparation method and analysis with fast GC-MS/MS could aid the forensic medicine and toxicology communities in elucidating important information from these often-overlooked biological matrices.

Biograph :

Brian Waters has received his Master of Science degree in Criminalistics from California State University, Los Angeles, USA. After working as a Criminalist for the County of Los Angeles, Department of Coroner/Medical Examiner for almost 8 years, he

has joined as an Assistant Professor in the Department of Forensic Medicine at Fukuoka University in Japan. His specialty is postmortem forensic toxicology and he has published academic papers on fast gas chromatography-mass spectrometry, the analysis of novel psychoactive compounds, preparation methods for postmortem samples and the analysis of volatile hydrocarbons in blood.